Bond Apatite®

Bone Graft Cement

Instructions for use

(EN) English

(FR) Français

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(IT) Italiano

(DE) Deutsch

(ES) Español

(RO) Română

(PT) Português

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Augma Biomaterials Ltd.

Alon Hatavor 20 St., P.O.Box 3089

Caesarea Southern Industrial Park

3088900, Israel

www.augmabio.com

Superseded

Bond Apatite®

Bone Graft Cement

Augma Biomaterials Ltd.
Alon Hatavor 20 St., P.O.Box 3089
Caesarea Southern Industrial Park
3088900, Israel
www.augmabio.com

English

INSTRUCTIONS FOR USE

Bond Apatite® is a bone graft cement composed of biphasic Calcium Sulfate and Hydroxyapatite. Please read this entire circular before performing the procedure. This device is for sale by, or on the order of a physician, or licensed practitioner.

COMPOSITION

Each driver contains highly pure medical grade biphasic calcium sulfate, hydroxyapatite granules and physiological saline.

PRODUCT DESCRIPTION

Bond Apatite® is a synthetic, osteoconductive, composite graft cement composed of biphasic calcium sulfate and hydroxyapatite granules. Bond Apatite® functions as a scaffold for bone regeneration in dental procedures, and is intended to fill, augment or reconstruct bony defects in the oral and maxillofacial region.

STERILIZATION

Bond Apatite® sterilization is carried out using γ - irradiation.

INDICATIONS FOR USE

Bond Apatite® is a synthetic osteoconductive, bone grafting cement composed of hydroxyapatite and biphasic calcium sulfate in granulated powder form, intended to fill, augment or reconstruct periodontal or bony defects of the oral and maxillofacial region. Examples of clinical applications:

- Sinus floor augmentation
- Periodontal bone defects
- Debridement; fenestrations
- Alveolar ridge augmentation
- Horizontal defect (and crest widening) filling bony defects pre implant placement
- Filling of cyst cavities

CONTRAINDICATIONS

• The customary contraindications in oral and maxillofacial surgery with other implant materials should be observed.

• Lack of adequate training of the practitioner is a major risk for the success of the implant procedure.

ABSOLUTE CONTRAINDICATION

Acute and chronic active infection at the site of the implant.

RELATIVE CONTRAINDICATION

• Serious disturbance of bone metabolism • Serious bone diseases of endocrine etiology • Severe or difficult to control diabetes mellitus • Severe renal dysfunction, severe liver disease • Vascular impairment at the implant site • Immunosuppressive and radiation therapy • Ongoing treatment with glucocorticoids, mineralocorticoids and growth agents affecting calcium metabolism • Malignancies • Lactating and pregnancy • Effect on pediatric patients is not known

PRECAUTIONS AND WARNINGS

• The material is supplied for one time use. DO NOT RESTERILIZE in case of partially opened or defective primary packing (driver) or secondary packing (peel-off blister). Under these circumstances Bond Apatite® may no longer be used since the sterility of the material is no longer ensured • Do not use when the temperature of the product is below 10°C (50°F). If used in low temperatures, wait until the product reaches room temperature (low temperatures will slow the setting reaction of the material) • It is not recommended to mix Bond Apatite® with other bone graft materials • It is not recommended to mix Bond Apatite® with blood • The expiration date is printed on the peel-off blister and on the packaging • Do not use after indicated expiry date • Bond Apatite® does not possess sufficient mechanical strength for load bearing defects prior to tissue in-growth. In cases where load support is required, standard internal or external stabilization techniques must be followed to obtain rigid stabilization in all planes • Bond Apatite® must not be used to stabilize screw placement • Bond Apatite® is intended for use by clinicians familiar with bone grafting procedures • Possible complications are the same as to be expected of autogenous bone grafting procedures. These may include: superficial wound infection, deep wound infection, deep wound infection with osteomyelitis, delayed union, loss of reduction, failure of fusion, loss of bone graft, graft protrusion and/or dislodgement, and general complications that may arise from anesthesia and/or surgery • Complications specific to oral/dental use are those as may be typically observed for similar bone grafting procedures and may include: tooth sensitivity, gingival recession, flap sloughing, resorption or ankylosis of the treated root, abscess formation, nerve injury.

STORAGE

Store at temperatures between 5°C (41°F) to 30°C (86°F). Avoid contact with a source of heat. Do not store the product in direct sunlight.

HANDLING

Bond Apatite® is available as granulated powder packed within a driver. The driver is separated into 2 compartments, one contains sterile standard saline liquid (0.9% Sodium Chloride for injection) and the second contains the powder.

IMPORTANT

• Read all steps of the instructions (1-3) before using Bond Apatite® • We highly recommend practicing the use of Bond Apatite® before first usage.
At the beginning of the procedure:

- Reflect the mucoperiosteal flap
- Remove the undesired soft tissue from the exposed bone surface.

Prepare the defected area for augmentation procedure
It is advisable to cover Bond Apatite® with a membrane barrier for optimal tissue regeneration.

STEP BY STEP


Read all instruction steps (1-3) before using Bond Apatite®. It is recommended that the time elapsed from the introduction of saline into Bond Apatite® should not exceed 3 minutes; otherwise the material will lose its pliability.

PREPARATION STEP

Before applying Bond Apatite®, make sure you have the items described:

A - A sterile Bond Apatite® driver

B - Dry sterile gauze pads

In order to achieve complete saturation of the powder, it is essential to flick the driver's cap gently a few times.

STEP 1.

Hold the driver's cap firmly and push the shaft slowly towards the line marked on the driver until you feel maximal resistance (this will activate the material and prepare it for ejection).

Note: While pushing the shaft, mild pressure is required.

STEP 2.

Remove the cap by twisting and pulling it out.

STEP 3.

Continue to push the shaft and inject the material into the required site.

Note: The paste should be in direct contact with the bone and slightly overfilled.

After ejecting the material into the required site:

Place a dry gauze pad on the material and condense the material firmly for 3 seconds.

Remove the gauze and slightly shape, if required.

Press on the material once again with dry gauze.

Remove the gauze and proceed with soft tissue coverage and wound closure.

POST TREATMENT CARE

- A surgical dressing may be placed over the wound for one to two weeks.
- Instructions that include an appropriate analgesic, antibiotic and home care regimen should be delivered to the patient after surgical intervention.
- When Bond Apatite® is used for filling bony defects pre implant placement at the site held for 3 months to 6 months before dental implant placement.
- Prior to implant placement the grafted site must be evaluated to ensure that adequate bone healing has occurred.
- Place the dental implant and abutment according to the clear indications and instructions for the dental implant/abutment system used.

RESORPTION PROFILE

The resorption of Bond Apatite® follows the resorption rate of the two components. The biphasic calcium sulfate component is completely resorbed in approximately 12 weeks.

The hydroxyapatite component is slowly resorbed in the body and may remain at the grafted site for months to years.

Do Not Use if Package is Damaged

Sterilized using gamma irradiation

Do not re-use

Use by date

Do not re-sterilize

Temperature sensitive

Date of manufacture (including sterilization)

Attention. See instructions for use

Catalogue number

Batch code

Consult electronic instructions for use

Caution:

US federal law restricts this product for sale by or on the order of a dentist or physician.

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